


Examination Couch

Constructed from chrome-plated square steel tube. Fully adjustable headrest. Top is upholstered and covered with washable plastic and flame retardant material

180 x 55 x 80 (H) cm
SOS 01-001



Gynecological Couch

Constructed from chrome-plated rectangular steel tube. Adjustable back rest and leg section.

Top is upholstered and covered with washable plastic and flame retardant material.

Complete with S.S. Steel bowl 32 cm and two adjustable knee crutches.

180 x 55 x 80 (H) cm
SOS 01-002



Gynecological Couch

Made of epoxy coated steel. Two section, one moveable have castors two of which have brakes and one fixed Complete with liquid S/S bowl holders.

135 x 92 x 86 (H) cm
SOS 01-003
68 x 92 x 86 (H) cm
SOS 01-004



Bed Panel Screen

Chrome-plated steel tubular frame, plastic curtains washable and frame retardent.

Mounted on castors.

2 Panels 100 x 42 x 167 (H) cm
SOS 01-005
3 Panels 150 x 42 x 167 (H) cm
SOS 01-006
4 Panels 200 x 42 x 167 (H) cm
SOS 01-007



HOSPITAL FURNITURE


Spygmomanometer Aneroid
on Castors



SOS 01-009




I.V. Stand on Castors
160 to 220 cm
SOS 01-008



I.V. Stand on Castors
160 to 220 cm
SOS 01-010

S.S Bottle Holder
SOS 01-013



003



S.S Four Hooks
SOS 01-011



S.S Two Hooks
SOS 01-012



Plastic Four Hook
SOS 01-014

HOSPITAL FURNITURE



Examination Lamp
SOS 01-015



Examination Lamp
SOS 01-016



Examination Lamp
SOS 01-017

Examination Lamp complete with Halogen bulb. Chrome Plated steel tubing. Foldable stainless steel base with 5 castors

X- Ray Film Viewers

Wall model. Constructed from enamelled sheet steel.

Flourescent cold light.

Unbreakable with monochromatic diffuser with supporting device.

40 x 43 cm	80 x 43 cm	120 x 43 cm
SOS 01-018	SOS 01-019	SOS 01-020



HOSPITAL FURNITURE

Instruments Trolley

Constructed from stainless steel 18/ 10. Shelves with upstands on four sides. Push handle. Mounted on 80 mm. Castors

60 x 40 x 80 (H) cm 70 x 50 x 80 (H) cm

SOS 01-021 SOS 01-022

90 x 60 x 80 (H) cm

SOS 01-023



Instruments Trolley

Constructed from stainless steel 18/ 10. Shelves with upstands on four sides. Push handle. Mounted on 80 mm. Castors

70 x 50 x 80 (H) cm

SOS 01-024



Trolley

Complete with S.S Bowl 18/10. 32 cm

70 x 50 x 80 (H) cm

SOS 01-025



Dressing Trolley

Complete with Stainless steel bottle rack

Rectangular waste bin S.S Bowl 32 cm

90 x 60 x 80 (H) cm

SOS 01-026



Bowl's Stand (Double)

Constructed from chrome plated steel tube
Variable height. Stainless steel bowl
Stainless steel Base (Double) SOS 01-027
Stainless steel Base (Single) SOS 01-028



Bowl's Stand (Single)

**Soiled Linen Trolley
Doubles**
SOS 01-030



Chrome plated frame.
Mounted on 80mm castors.
Removable bag.

Mayo Table

Stainless Steel 18/ 10 construction,
variable height from 80 to 115 cm
4 castor base 50 mm
Removable stainless steel tray
70 x 45 x 3 (H) cm
SOS 01-029



**Soiled Linen Trolley
Singles**
SOS 01-031





HOSPITAL FURNITURE



Mounting Step

Chrome plated steel frame with plastic feet

Step covered with anti slip rubber

1 Step SOS 01-032

2 Steps SOS 01-033





Kick Bucket

Manufactured from 18/ 10 stainless steel. Mounted on 5 castors 40 mm

41 x 18 H


SOS 01-034



SOS 01-039

Antitrauma Scooping Stretcher

SOS 01-040




Foldable in 2 (Length)

Open size 203 x 50 x 14

SOS 01-035

Closed size 96 x 50 x 14

SOS 01-036



Foldable in 4 (Length)

Open size 203 x 50 x 15

SOS 01-037

Closed size 96 x 17 x 15

SOS 01-038

Antitrauma Scooping Stretcher

Minimum length	167cm
Maximum length	210cm
Folded length	120cm
Width	43cm
Height	9cm
Weight	10kg
Load	170kg